



ECG BD 351 Foam Soap dispenser

Product code: 100000969041

EAN: 8592131307971

- Touch-free operation minimizing the spread of bacteria
- 2 different dose volume settings
- "Anti Drip" function
- IPX 4 waterproof

- Ideal for use in the kitchen and bathroom
- Touch-free operation minimizes the spread of bacteria
- "Anti Drip" function – prevents spontaneous dripping
- Waterproof IPX 4
- Automatic dosing system
- Precise scanning with infrared sensor

- Easy refilling
- 2 different dose volume settings (0.7 g / 1.4 g)
- Nozzle angle for convenient and accurate dosing
- LED indicator:
 - Blue color Switch-on / Dose setting approx. 0.7 g
 - White color / Dose setting approx. 1.4 g
 - Red color Switch off / Low battery indicator

- Only suitable for viscous liquid soap (diluted in ratio): 1 (liquid soap) : 2–6 (water)
- White finish
- Wall mounting option (mounting kit included)
- Diluted soap dispenser volume: 350 ml
- Power: 3× 1,5 V AA batteries (not included)

Technical specifications

Power supply	Battery
Type of battery	3x1,5V AA

Product specification

Liquid Soap	Yes
Foam Soap	No
LED indicator light	Yes
Sensor range	< 4 cm
Volume (ml)	350
Liquid amount adjust	2
Drop volume in (ml)	0,7 / 1,4

No Drip System	Yes
Water resist level	IPX4
Ergonomically shaped	Yes

Design

Colour	White
Material	Plastic (ABS)

Accessories

Drip tray	No
Protective cover	No
Battery included	No

Weight & dimensions

WIDTH (MM)	70
HEIGHT (MM)	205
LENGTH (MM)	85
WEIGHT (KG)	0.3

Weight & dimensions (Brutto)

WIDTH WITH PACKAGING (MM)	103
HEIGHT WITH PACKAGING (MM)	232
DEPTH WITH PACKAGING (MM)	84
WEIGHT WITH PACKAGING (KG)	0.8

Minimize bacterial spread

The **ECG BD 351 Foam** non-contact soap dispenser fits perfectly into any modern household, whether in the bathroom or the kitchen. The non-contact use of the dispenser minimizes the spread of bacteria, which will be appreciated especially by families with children. The infrared sensor takes care of starting the dispenser at the right moment and the automatic dosing system perfectly measures the dose of liquid soap. Depending on your preferences, you can choose from two dose sizes, 0.7 g or 1.4 g.



The device has an extra large container and can hold up to **350 ml** of liquid soap. So it will last you long enough and its topping up is very simple. We recommend using viscous liquid soap diluted in a ratio of 1 (soap): 2–6 (water).

Colour information indicators

The ECG BD 351 Foam dispenser is also equipped with practical **light indicators**, so that you will always be informed about the current setting. The blue light indicates that the dose is switched on and at the same time its setting to approx. 0.7 g. White means that you have set the dose to 1.4 g. Red indicates that it is switched off and also warns you of a low battery.



It does not drip and is waterproof

You will appreciate the popular **Anti Drip** function, which prevents unplanned dripping. You will get rid of unwanted smudges on the sink left by other similar dispensing systems. Also, the advantage of the ECG dispenser is its **waterproofing** IPX 4, so that you can be sure that it withstands water splashing.



The package includes a practical set for wall mounting. So you can choose where to place the dispenser. It is powered by 3× 1,5V AA batteries (not included).